

The background features a complex, abstract design with overlapping geometric shapes in various shades of blue (from light cyan to deep navy) and black. The shapes include triangles, circles, and curved segments, creating a dynamic and modern visual effect.

ISOCOAT

UV BLOCKER XL

User Manual

ISOCOAT

UV BLOCKER XL

User Manual

PRODUCT DESCRIPTION

ISOCOAT UV BLOCKER XL is a universal, light-curing gel resistant to high temperatures, designed for precise laboratory work in prosthetics and orthodontics — suitable for both plaster and 3D-printed models.

BEFORE USE

- This manual provides information on proper use and safety precautions to ensure safe and effective application of the product.
- Before the first use, allow the gel to reach room temperature.
- Do not expose the syringe to sunlight or UV radiation. Store the product at room temperature, away from heat sources and open flames.
- Make sure the surface is properly cleaned and cured. Apply the gel only when it is dry and dust-free.
- Do not use intraorally.

APPLICATION

- Apply the desired amount of gel directly from the syringe onto the area to be blocked. You may spread it using the supplied applicator, a brush, or any laboratory instrument.
- Remove any excess gel. For model repair, you can adjust the shape with a bur after polymerization.
- Cure the gel under any UV light source until fully hardened. It is compatible with most types of lamps — including handheld dental curing lights, UV print curing lamps, and nail curing units.
- After complete polymerization, wipe off the inhibition layer using a paper towel.

HAZARD IDENTIFICATION

- ISOCOAT UV BLOCKER XL may cause irritation. Use appropriate personal protective equipment.
- For detailed hazard and safety information, refer to the product label and Safety Data Sheet (SDS).

LOT NUMBER/ EXPIRY DATE

- The product's expiration date is printed on each ISOCOAT package.
- When contacting technical support, always provide the product's batch number.
- Do not use the product after its expiration date.


LABEL SYMBOLS



See the information about the insulator on the website www.syntegre.com



Store in a dark, dry place away from sun and heat sources.



Store the gel at the right temperature.

CONTACT INFORMATION

For product support, refer to the information available to you on the website: www.syntegre.com or write to us at: support@syntegre.com



Syntegre s.c.
ul. Międzyrzecka 41/36
21-400 Łuków
www.syntegre.com info@syntegre.com

rev. 2024-04