

ISOCOAT

UV BLOCKER XL

Technical Data Sheet

PRODUCT DESCRIPTION

ISOCOAT UV BLOCKER XL It is a universal, lightcuring gel resistant to high temperatures, designed for precise laboratory work in prosthetics and orthodontics — suitable for both plaster and 3D-printed models.

It works perfectly for modeling spaces in whitening trays, blocking undercuts and other detailed laboratory cases. Thanks to its optimized viscosity and thixotropy, the gel stays exactly where it's applied — it doesn't flow or spread, giving full control over the thickness and shape of the layer. Its transparent blue color ensures excellent visibility during application. High temperature resistance provides stability not only during acrylic polymerization and thermoforming, but also during the injection of thermoplastic materials.

AVAILABLE PACKAGES AND COLOURS

PRODUCT	COLOUR	SIZE	INDEX
ISOCOAT UV BLOCKER XL	SKY BLUE	3ml / 4g	IC1003
ISOCOAT UV BLOCKER XL	SKY BLUE	4 x 3ml / 4g	IC1004

Our products are continuously developed and researched, and we reserve the right to change product parameters without prior notice.